



[05] INTENSIVE AGE DEFENSE

Environ's Youth EssentiA Range is an antiaging skincare regimen that combines vitamins, antioxidants and peptides to create healthier, younger-looking skin. The line includes: Hydra-Intense Cleansing Lotion that removes makeup and impurities; moisturizing and antiaging Vita-Peptide Toner; fast-absorbing Vita-Peptide Eye Gel, which improves skin tone and elasticity; Antioxidant Defense Crème to hydrate and fight free radicals; and Vita Peptide C-Quence Serum, available in four levels of increasing concentrations of vitamins A, C and E.

Contact: 508.539.8900, dermaconcepts.com.

[06] ADD-ON SUPPORT

The PCA SKIN Advanced Treatment Booster is a professional retinoid treatment that may be used on its own or in conjunction with other peels or treatments. Key ingredients include: a vitamin A (10.5%) blend that promotes a clear complexion and evens skin tone; lentil seed extract (3%), which restores keratinization after treatment, limits pore wall slacking and reduces sebum production; soothing panthenol, escin, butcher's broom, gotu kola extract and marigold flower extract (2%); and vitex agnus-castus extract, tocopherol and ascorbyl tetraisopalmitate for hydration and antioxidant protection.

Contact: 877.722.7546, pcaskin.com.



[07] SOFTENING OIL


Soft Touch Facial Oil from Skinprint's Eco-Ceuticals line is made of Tahitian Monoi oil infused with coconut, rosehip and linseed oils. The nourishing ingredients tone, hydrate, soothe, soften and renew sensitive and aging skin, and the lightweight oil immediately absorbs into the skin without a greasy feel.

Contact: 800.234.1308, skinprint.com.

[08] PRE- AND POST-CARE

NEOCUTIS PREVedem Bruise Minimizing Serum helps diminish bruising and puffiness, and strengthens skin before and after cosmetic procedures. Patent-pending Hamamelis Virginiana Extract provides antioxidant support to help prepare skin for stress caused by procedures, a blend of polyphenols and essential nutrients nourish and prime skin to resist the appearance of bruising, and glycerin hydrates and soothes. Recommended use is twice daily for two weeks prior to aesthetic treatments—including injections, lasers, dermabrasion and peels—and continued application for one week afterward.

Contact: 866.636.2884, neocutis.com.



[09] EASIER ASPIRATION

The Lee Aspirator is a patent-pending safety device that snaps onto cosmetic injectable syringes, allowing providers to aspirate the syringe with one hand during dermal filler injections. The customized Thumb Ring lets you pull back on the plunger with minimal movement, and the Finger Clip prevents your fingers from sliding during aspiration. The reusable aspirator fits all Juvéderm syringes and is being designed to fit the new Restylane syringe.

Contact: 800.231.7762, air-tite.com.